

Keops User Guide

1st Edition

April 2007



Keops User Guide

Copyright © 2007 Matisse Software Inc. All Rights Reserved.

Matisse Software Inc.
930 San Marcos Circle
Mountain View, CA 94043
USA

Printed in USA.

This manual and the software described in it are copyrighted. Under the copyright laws, this manual or the software may not be copied, in whole or in part, without prior written consent of Matisse Software Inc. This manual and the software described in it are provided under the terms of a license between Matisse Software Inc. and the recipient, and their use is subject to the terms of that license.

RESTRICTED RIGHTS LEGEND: Use, duplication, or disclosure by the government is subject to restrictions as set forth in subparagraph (c)(1)(ii) of the Rights in Technical Data and Computer Software clause at DFARS 252.227-7013 and FAR 52.227-19.

The product described in this manual may be protected by one or more U.S. and international patents.

TRADEMARKS: MATISSE and the MATISSE logo are registered trademarks of Matisse Software Inc. All other trademarks belong to their respective owners.

Content

INTRODUCTION.....	8
SCOPE	8
BEFORE YOU BEGIN.....	8
USING KEOPS.....	8
GETTING STARTED.....	10
CONFIGURING KEOPS.....	10
LAUNCHING KEOPS.....	10
CONFIGURING KEOPS.....	11
CONFIGURING THE DATABASE SERVER.....	11
CONFIGURING A WORKSTATION	11
<i>Setting up the Database Server Connection</i>	11
<i>Defining an Administrator</i>	11
<i>Defining Users</i>	12
<i>Setting up Lexicons</i>	12
<i>Configuring a Consultation Workstation</i>	12
<i>Configuring an Acquisition Workstation</i>	12
MANAGING USERS.....	13
USER ACCOUNT TYPES	13
USER ROLES	13
VIEWING USERS.....	13
ADDING A USER.....	13
UPDATING A USER	14
DELETING A USER.....	14
MANAGING MEDIA ACQUISITION INSTRUMENTS.....	15
MEDIA DEVICE MODELS.....	15
<i>Canon EOS Digital Camera</i>	15
<i>FlashBus Spectrim Video Board</i>	15
<i>Sony 1394 Digital Video Camera</i>	15
VIEWING INSTRUMENT SETUPS.....	15
SETTING UP AN INSTRUMENT	16
UPDATING AN INSTRUMENT SETUP	16
DELETING AN INSTRUMENT SETUP	16
CHECKING MEDIA CAPTURE	16
<i>Digital Still Camera-based Instrument</i>	17
<i>Video Camera-based Instrument</i>	17
CREATING PREDEFINED VIDEO SETTINGS	17
UPDATING PREDEFINED VIDEO SETTINGS.....	19
DELETING PREDEFINED VIDEO SETTINGS	19

MANAGING USER-DEFINED LEXICONS.....	20
LEXICON-BASED PROPERTIES	20
SETTING UP A LEXICON.....	21
MANAGING STUDIES	22
VIEWING STUDIES.....	22
CREATING A NEW STUDY	22
UPDATING A STUDY.....	22
DELETING A STUDY	22
MANAGING OBSERVED ANIMALS	24
VIEWING ANIMALS	24
CREATING AN NEW ANIMAL.....	24
UPDATING AN ANIMAL	24
DELETING AN ANIMAL.....	24
MANAGING OBSERVATIONS.....	26
VIEWING OBSERVATIONS.....	26
RECORDING AN OBSERVATION	26
<i>Managing Animals</i>	26
<i>Managing Imaging Instruments</i>	27
<i>Choosing Pictures</i>	28
<i>Creating an Observation</i>	28
<i>Deleting Captured Images</i>	29
UPDATING AN OBSERVATION.....	29
DELETING AN OBSERVATION	29
COMPARING OBSERVATIONS	29
VIEWING IMAGE DETAILS	30
EXPORTING AN IMAGE	30
QUERYING THE STUDIES DATABASE	31
BUILDING A QUERY	31
EXPORTING A QUERY RESULT	31
MANAGING THE AUDIT RECORDS.....	33
RECORDED OPERATIONS.....	33
VIEWING THE AUDIT LOG	33
EXPORTING THE AUDIT LOG	33
SESSION EVENT LOG.....	35
VIEWING SESSION EVENTS	35
ACCESSING SESSION EVENTS HISTORY	35
MANAGING THE DATABASE SERVER.....	36
CREATING THE KEOPS DATABASE	36
<i>Configuring the Database</i>	36

<i>Starting the Database</i>	36
<i>Creating the Database User</i>	36
<i>Loading the Database Schema</i>	37
STOPPING THE DATABASE SERVER	37
RESTARTING THE DATABASE SERVER.....	37
BACKING UP THE DATABASE	37
EXPORTING AND IMPORTING DATA.....	37
<i>Exporting Data</i>	37
<i>Importing Data</i>	38
MIGRATING STUDIES.....	40
DATA MODEL	40
FILE FORMAT	40
EXPORTING DATA	41
LOADING DATA	42
<i>Re-initialize the Database</i>	42
<i>Load Keops Database Schema</i>	42
<i>Load Keops Migration Schema Extensions</i>	42
<i>Load CSV Files</i>	42
FREQUENTLY ASKED QUESTIONS.....	44
CONFIGURATION	44
USERS	44
ACQUISITION INSTRUMENTS	44
USER-DEFINED LEXICONS	44
STUDIES	45
AUDIT RECORDS	45

Table of Figures

FIGURE 1: SUPPORTED VIDEO FEATURES	18
FIGURE 2: SONY 1394 CONTROL FEATURES	18
FIGURE 2: PROPERTIES MANAGED BY A USER-DEFINED LEXICON	20
FIGURE 3: EXAMPLE OF OBSERVATION SELECTIONS PER CATEGORY VALUE	21
FIGURE 4: PROPERTY TYPE FORMATS	31
FIGURE 5: AUDITED OPERATIONS	33
FIGURE 6: STUDY DATA MODEL	40

Introduction

Keops provides a secure, efficient and easy to use solution to collect, analyze, process and store high-quality digital media of pre-clinical or clinical studies.

Keops features a distributed system to manage digital assets of pre-clinical or clinical studies. The main features of Keops are comprised of:

- Connectivity to imaging instruments to capture digital media. Digital still cameras, analog video camera and digital video cameras are currently supported.
- Management of pre-clinical or clinical studies with several treatment groups.
- Management of users and modifications in compliance with FDA regulations. All changes related to clinical study data are recorded and can be viewed in an audit log.
- A comprehensive set of in-depth searching capabilities.

Scope

This document is intended to help new users learn the configuration, administration and utilization of Keops features.

Before you begin

Throughout this document, we presume that you already know how to manipulate the external imaging instruments connected to Keops. We also assume that Keops installation has been successfully completed.

Using Keops

Keops User Interface is designed to allow users to quickly execute repetitive operations using only the workstation keyboard.

The main dialog windows appear by pressing a keystroke combination on the keyboard. For example, pressing `Ctrl+Shift+S` on the keyboard displays the *Add Study* dialog window.

Each capture form dialog window is designed to quickly complete the form by moving the cursor from one field to the next in the most convenient order for

the operator. Pressing `Tab` on the keyboard moves the cursor from fields to fields up to the OK button. Pressing `Shift+Tab` on the keyboard moves the cursor backwards.

When a menu selection is activated, pressing the first letter of a selection on the keyboard selects it. The `Up` and `Down` arrows on the keyboard move selections respectively up and down. Pressing `F4` displays all possible selections.

When an image thumbnail list viewer is activated, use the arrows (`Up`, `Down`, `Right` and `Left`) on the keyboard to move from one image thumbnail to the next. When the list viewer enables multi-selection, use `Shift+Arrow` to select multiple image thumbnails.

When a button is activated, pressing `Enter` on the keyboard validates the selection.

When an image, a thumbnail or a row in a query result table is selected, right-click on the mouse to get access to a contextual menu.

Getting Started

Configuring Keops

If you have just completed the installation, you must start from the *Configuring Keops* section of this document.

Launching Keops

To launch Keops, select `Keops` in the `Matisse Life Science` folder in the `Start All Programs` menu on your machine.

NOTE: You will be prompted for login and password to enter into the application.

Configuring Keops

Keops is a distributed system that requires the setup of a database server and of one or more workstations.

CAUTION: You cannot configure a workstation unless the database server is configured.

Configuring the Database Server

The *Creating the Keops Database* section of this document guides you through the database configuration process. You need to refer to the *Creating the Keops Database* section and to complete the instructions presented in this section.

Once you have successfully completed the database server setup, you can move on to configuring a workstation.

Configuring a Workstation

Setting up the Database Server Connection

The very first time you launch Keops on a workstation, the *Database Connection Setup* dialog window appears so you can provide the necessary information to establish a connection to the database server.

The `Ping` button helps you check the database connection prior to saving the database connection settings.

NOTE: This operation is done only once per workstation since all users on a workstation are sharing the database connection settings.

Defining an Administrator

The very first user who logs into the application becomes an administrator. When you have successfully passed the *Login* dialog window, the *Create a User* dialog window is displayed and the administrator full name can be set. You are now logged in as an *administrator* ready to complete the setup.

NOTE: From now on, this administrator is the only user able to login into Keops until you have defined new users.

Defining Users

Since you are now logged in as an *administrator*, you can manage application users. The *Managing Users* section of this document describes the management of users. You can add new users now or later.

Setting up Lexicons

Since you are now logged in as an *administrator*, you can manage user-defined application lexicons. Section *Managing User-defined Lexicons* of this document describes the management of application lexicons. You can configure the lexicons now or later, but since Keops is not pre-configured with default selections, you need to define domain-related selections for each property prior to using Keops in production.

Configuring a Consultation Workstation

Once you have completed the database connection setup, you need to connect as an *administrator* and check the data is accessible before quitting the application. This completes the consultation workstation configuration.

Configuring an Acquisition Workstation

Once you have completed the database connection setup, you need to connect as an *administrator* to configure the acquisition instruments connected to this workstation.

The *Managing Media Acquisition Instruments* section of this document guides you through configuring acquisition instruments.

Managing Users

To manage the users, log in as an *administrator* and select the `Users` sub-tab under the `Administration Tab`.

User Account Types

The application manages 2 types of user account:

System User: a username registered into the application, which matches the Windows login name of the operator. No password is registered into Keops since the application checks against Windows credentials.

Application User: a username and a password registered into Keops. The username does not need to match any Windows login name.

User Roles

The application manages 2 types of role:

Operator: The operator main role consists in managing studies. Specific tasks include creating studies, animals, capturing photos and establishing diagnostics associated to the captured photos. The operator can also initiate the change or removal of existing study information, but an approval by an administrator is required to complete the operation.

Administrator: the administrator can create, update and delete *System* and *Application* users. He can also grant and revoke privileges to the users. He can also set and reset the password of any application users.

The administrator can also configure the application. He manages instrument setups as well as application lexicons.

The administrator is also responsible for approving changes or removals initiated by an operator.

Viewing Users

The application users are listed in the `Users` sub-tab under the `Administration Tab`.

Click on the `Select a Username` scrolling list to view the detailed information about an application user.

Adding a User

To add a new user, click on the `Create new user` button in the toolbar. The *Create a User* dialog window appears. You can then enter a new username and

define the user's full name, role and connection mode and enter a note. The username value must be unique. Optional fields include passwords and note.

An entry in the audit log is added to record the operation.

NOTE: When a new user is added to Keops, it becomes immediately visible to all Keops workstations.

Updating a User

To update a user, click on the `Update selected user` button in the toolbar. The *Update a User* dialog window appears so you can change the user full name, role and connection mode.

An entry in the audit log is added to record the operation.

Deleting a User

To delete a user, click on the `Delete selected user` button in the toolbar. The *Delete a User* dialog window appears allowing you to confirm your selection before proceeding with the removal of the user.

An entry in the audit log is added to record the operation.

Managing Media Acquisition Instruments

Acquisition Workstations are dedicated to the capture of digital media assets. These workstations are directly connected to the imaging instruments that capture digital media.

To manage the instrument setups, log in as an *administrator* and select the `Instruments` sub-tab under the `Administration Tab`.

Media Device Models

Keops supports connectivity to the following 3 types of imaging devices:

Canon EOS Digital Camera

Keops supports the Canon EOS Digital SLRs with DIGIC II processors connected for example to a non-mydratiac digital fundus camera (e.g. a Canon EOS20D connected to a Canon CR66Dgi).

FlashBus Spectrim Video Board

Keops supports the FlashBus Spectrim product family, a high-performance, PCI bus-mastering video frame grabber designed to capture composite, S-Video, or RGB video in real time. For example, the video board can be connected to an analog video indirect ophthalmoscope.

Sony 1394 Digital Video Camera

Keops supports the Sony “FireWire” digital video cameras (Color models only) including the XCD-SX910CR/X710CR series. For example, a Sony 1394 digital video camera can be mounted on a digital video slit lamp.

An acquisition workstation can be connected to at most 3 instruments. Only one instrument for each media device model can be connected to a workstation at a time.

Viewing Instrument Setups

The connected instruments are listed in the `Instruments` sub-tab under the `Administration Tab`.

Click on the `Select an InstrumentID` scrolling list to view the detailed information about an instrument.

Setting up an Instrument

CAUTION: The external device drivers for the instrument must be installed prior to complete the instrument connection setup.

To setup a new instrument connection, click on the `Setup new instrument` button in the toolbar. The *Setup a new Instrument* dialog window appears. You can then select the media device model and enter a new instrument Id, serial number, location and note. The Instrument Id value must be unique. Optional fields include serial number, location and note.

Additional options are provided in the `Advanced` Tab when you setup the connection to a Sony 1394 “FireWire” digital video camera. You can select a video format ranging for high image resolution (1280x960) to low image resolution (640x480). By default, the video format is set to the highest resolution. You can also configure a serial port for switch-based capture of digital images. An action on the external switch (i.e. a foot pedal or a joystick) will have the same effect as clicking on the `Take a Picture` button in the toolbar. By default, the switch-based capture is not enabled.

An entry in the audit log is added to record the operation.

Updating an Instrument Setup

To update an instrument setup, click on the `Update selected instrument` button in the toolbar. The *Update Instrument* dialog window appears allowing you to update the instrument serial number, location and notes. The properties defined in the `Advanced` Tab can also be updated.

An entry in the audit log is added to record the operation.

Deleting an Instrument Setup

To delete an instrument setup, click on the `Delete selected instrument` button in the toolbar. The *Delete an Instrument* dialog window appears so you can confirm your selection before proceeding with the removal of the instrument setup.

An entry in the audit log is added to record the operation.

Checking Media Capture

Once you have checked that the equipment is correctly plugged to the workstation and the equipment power is turned on, click on the `Connect to selected instrument` button from the toolbar, the status in the *Live Feed* section of the window is now indicating `Connected`. If not, you need to check the event session log for more information about the connection issue. If an external switch has been configured the status also indicates `External switch active`.

The instrument is now connected and some buttons in the toolbar are enabled. There are 2 scenarios depending on the media device type mounted on the instrument.

Digital Still Camera-based Instrument

The instrument uses a Canon EOS Digital camera, which is now connected to Keops. The `Disconnect` button in the toolbar is now enabled while the `Connect` button is disabled. The `Take a Picture` button is also enabled. The `Download images` button may be enabled as well if there are images stored in the camera waiting to be downloaded. Images are stored in the camera when the camera is on and that pictures are taken while the instrument is not connected to Keops.

Pictures can be taken by clicking on the `Take a Picture` button or by pressing the shutter button of the digital camera or an external switch plugged into the camera.

Video Camera-based Instrument

The instrument uses a video camera, which is now connected to Keops. The `Disconnect` button in the toolbar is now enabled while the `Connect` button is disabled. The `Start Video` button is also enabled. Click on the `Start Video` button to capture the live video feed.

When you see the live video feed in the *Live Feed* section of the window, you can take pictures by clicking on the `Take a Picture` button or by pressing the external switch if the status indicates that one switch is active.

When the video is live, you can also define new video settings by clicking on the `Create, Update or Delete video settings` buttons in the toolbar.

Creating Predefined Video Settings

CAUTION: Operations on video settings are available only for video-based instruments and are only enabled when the video is live.

The objective of the predefined video settings is to speed up the image capture during the observations. For the same instrument, you can define a set of video settings well suited to each capture situation.

Keops manages all video features supported by the imaging device. The table below summarizes the features supported by each video-based device.

Features	Sony Video Camera	FlashBus Video Board
Brightness	Yes	Yes
Contrast	No	Yes
Saturation	No	Yes
Hue	Yes	Yes
Shutter	Yes	No
Gain	Yes	Yes
White Balance	Yes	No

Figure 1: Supported Video Features

The control features of the Sony “FireWire” digital video camera details as follows:

Features	
Brightness	Adjust the black level precisely
Shutter	Controls the exposure time. Can be controlled by both: Relative control values from 1/100,000 of a second to 16 seconds allocated from 3 to 1150. Setting examples: 3: 10 micro sec 32: 1 ms 100: 10 ms 1000: 1 s 1010: 2 s 1050: 16 s
Gain	Absolute value control from values 1/100,000 of a second to 17.5 seconds set continuously.
White Balance	Can be changed to 0 to 18 dB, subdivide in 640 steps
Hue	Adjust the White balance by adjusting the R and B gain with respect to G Adjust G gain. Use this feature when you cannot obtain the correct white balance using the R and B gain

Figure 2: Sony 1394 Control Features

CAUTION: Long exposure times will slow down the frame rate in accordance with the exposure time. Shutter above 300 have a significant impact on the frame rate.

To setup new video settings, click on the `Create new video settings` button in the toolbar. The `Create new settings` dialog window appears. You can then enter the settings name, define the default picture shot count and adjust the video filters in real-time on the live feed.

Updating Predefined Video Settings

To update the video settings, click on the `Update selected video settings` button in the toolbar. The *Update Settings* dialog window appears allowing you to change the default picture shot count and to adjust the video filters in real time on the live feed.

Deleting Predefined Video Settings

To delete the video settings, click on the `Delete selected video settings` button in the toolbar. The *Delete Settings* dialog window appears so you can confirm your selection before proceeding with the removal of the video settings.

Managing User-defined Lexicons

Keops allows you to tailor property value selections to your domain of application. Keops is not pre-configured with default selections therefore you need to define domain-related selections for each property prior to using Keops in production.

Lexicon-based Properties

Keops enables to define lexicons of closed selections for the following study entity properties:

Entity	Property
Study	Species
Animal	Group Sex
Diagnostic	Laterality Category Observation

Figure 3: Properties managed by a user-defined lexicon

Properties are usually not linked of each other thus selections for one property do not affect the selections for another property. For example, a selection in property *Group* is not affecting the selections in property *Sex*. However, the *Category* and *Observation* properties can be linked allowing you to define a sub-lexicon of Observation selections associated to one Category selection. For example as shown in the table below, *Cornea* and *Iris* may use specific Observation selections while others such as *eyelid*, *lens*, *retina*, etc. are using the default selections.

Category	Observation
Cornea	focal opacity pinpoints apocities ulcer
Iris	vascularization œdème
All other selections	focal opacity pinpoints apocities ulcer vascularization œdème anterior synechia

Figure 4: Example of Observation selections per Category value

Setting up a Lexicon

To manage lexicons, log in as an *administrator* and select the `Lexicons` sub-tab under the `Administration` Tab.

Select an Entity and a Property and then add new selections to a list. The list of selections is displayed in the order you have organized them. Each lexicon must hold at least one entry besides the selection for “no value”.

NOTE: The “no value” selection is managed by the application therefore there is no need to define such a selection in the lexicon.

NOTE: The selections are text-based and the selected selections are recorded as is, but comparison between 2 selections is not case sensitive.

CAUTION: Keops is not pre-configured with default selections therefore you need to define selections for each property prior to using Keops in production.

CAUTION: An administrator can update property selections at any time, since an operator can only choose from the selections pre-defined in the lexicon.

Managing Studies

In Keops, a Study is defined by a set of observed animals for which a set of observations has been recorded. A picture shot is associated to each recorded observation.

To manage studies, select the `Studies` Tab.

Viewing Studies

The studies are listed in the `Studies` Tab.

All recorded studies are listed in the `Select a Study Number` menu selection. Click on the `Select a Study Number` scrolling list to view the detailed information about a study. Then to navigate through the hierarchy of data, click on the `Select an Animal ID` scrolling list to switch from one animal to the next in the list. Then for a selected animal, the observations are listed in an image list viewer. Selecting an examination presents the observation details including photos and diagnostics.

Creating a new Study

To create a new study, click on the `Create new study` button in the toolbar. The *Create a Study* dialog window appears. You can then enter a new study number, compound and select the species, a start date for the study and enter a note. The study number value must be unique. Optional fields include compound, species and note.

Updating a Study

To update a study, click on the `Update selected study` button in the toolbar. The *Update Study* dialog window appears allowing you to update the study number, compound, species, start date and notes. Each property update requires the approval of an administrator.

An entry in the audit log is added to record the change of each property.

Deleting a Study

To delete a study, click on the `Delete selected study` button in the toolbar. The *Delete a study* dialog window appears so you can confirm your selection before proceeding with the removal of the study. Deleting a study requires the approval of an administrator.

An entry in the audit log is added to record the operation.

CAUTION: When deleting a study all the depending entities are deleted as well. This includes the animals, exams, diagnostics and images part of this study.

Managing Observed Animals

In Keops, observed animals are uniquely referenced in a study. Each observed animal can carry a set of observations.

To manage observed animals, select the `Studies` Tab.

Viewing Animals

Observed animals, which are parts of studies are listed in the `Studies` Tab.

Once you have selected a study from the `Select a Study Number` scrolling list, click on the `Select an Animal ID` scrolling list to switch from one animal to the next. Then for a selected animal, the observations are listed in an image list viewer. Selecting an examination presents the observation details including photo and diagnostics.

Creating an new Animal

To create a new observed animal, click on the `Create new animal` button in the toolbar. The *Create an Animal* dialog window appears. You can then enter a new initial Id, select the sex and group and enter a study Id and a note. The animal initial Id value must be unique. Optional fields include sex, group, study Id and note.

Updating an Animal

To update an observed animal, click on the `Update selected animal` button in the toolbar. The *Update Animal* dialog window appears allowing you to update the animal initial Id, sex, group, study Id and note. Each property update requires the approval of an administrator.

An entry in the audit log is added to record the change of each property.

Deleting an Animal

To delete an observed animal, click on the `Delete selected animal` button in the toolbar. The *Delete an Animal* dialog window appears so you can confirm your selection before proceeding with the removal of the animal. Deleting an observed animal requires the approval of an administrator.

An entry in the audit log is added to record the operation.

CAUTION: When deleting an observed animal all the depending entities are deleted as well. This includes the exams, diagnostics and images part of this animal.

Managing Observations

In Keops, observations are uniquely referenced in an observed animal.

To manage observations, select the `Studies` Tab.

Viewing Observations

Observations that are parts of Observed animals are listed in the `Studies` Tab. Once you have selected an observed animal from the `Select an Animal ID` scrolling list, you can view the observations, which are listed in an image list viewer. Selecting an observation presents the observation details including photos and diagnostics.

Recording an Observation

An Observation recorded in Keops always contains a picture to support the operator diagnostic. The process for recording an animal observation implies that the operator takes meaningful pictures to document its findings.

To record an observation, click on the `Capture Media` button in the toolbar. The *Capture Media* dialog window appears. The operator can then capture images of an animal under observation.

From the *Capture Media* window, you can independently manage observed animals and imaging instruments.

Managing Animals

Observed animals, which are parts of studies are listed in the `Animal Initial ID` scrolling list. When the *Capture Media* window is opened, the selected animal is the one selected in the `Study` Tab. The operator can switch animals by clicking on the `Animal Initial ID` scrolling list.

If the animal is not yet present in the list, click on the `Create new animal` button from the toolbar to create a new animal.

In some cases, the operator may capture images of an animal that is not yet registered into the study and then decide after the picture shots examination that the animal needs to be added to the study. In this context, the operator needs to select the `Unregistered` value in the `Animal Initial ID` scrolling list or needs to click on the `Observe an unregistered animal` button from the toolbar.

If the `Unregistered` animal is selected, the operator will be prompted to create a new animal prior to recording an observation. And the remaining captured images will be associated to the newly created animal.

CAUTION: Captured images associated to an `Unregistered` animal are not shared with other workstations. These images are automatically deleted when the Keops session is ended.

Managing Imaging Instruments

Imaging instruments connected to the acquisition workstation are listed in the `Instrument` scrolling list. When the *Capture Media* window is opened, the selected instrument is the first one in the `Instrument` scrolling list. The operator can switch instruments by clicking on the `Instrument` scrolling list. Next to the `Instrument` scrolling list is displayed the status of the instrument.

Click on the `Connect` and `Go live` button from the toolbar, the instrument status is now indicating `Connected` or `Live`. If not, you need to check the event session log for more information about the connection issue. If an external switch has been configured, the status also indicates `External switch active`.

The instrument is now connected and some buttons in the toolbar are enabled. There are 2 scenarios depending on the media device type mounted on the instrument.

Digital Still Camera-based Instrument

The instrument uses a Canon EOS Digital camera, which is now connected to Keops. The `Take a Picture` button in the toolbar is enabled. The `Download images` button may be enabled as well if there are images stored in the camera waiting to be downloaded. Images are stored in the camera when the camera is on and that pictures are taken while the instrument is not connected to Keops.

Pictures can be taken by clicking on the `Take a Picture` button or by pressing the shutter button of the digital camera.

Video Camera-based Instrument

The instrument uses a video camera, which is now connected to Keops. The `Take a Picture` button in the toolbar is enabled. The `Take a Picture and Save` button is also enabled as well as the `Adjust video settings` button.

When you see the live video feed in the *Live Feed* section of the window, you can take pictures by clicking on the button or by pressing the external switch if the status indicates that one switch is active.

When the video is live, the operator can switch the video settings by clicking on the `Settings` scrolling list. Video settings can be adjusted in real-time by clicking on the `Adjust video settings` button in the toolbar.

NOTE: We recommend that an administrator check the instrument connectivity and settings prior to start the observations recording session.

Choosing Pictures

When multiple images of the same observation have been taken, the operator may need to review the images in detail and compare them before deciding for the more relevant picture shots.

To review and compare images, click on the `Select Media` button in the toolbar. The *Select Media* dialog window appears so the operator can review and compare images side by side.

To help the operator with the comparison, selected images can be displayed in single, dual or quad mode.

Creating an Observation

Saving a captured image displays the *Add an Exam* dialog window allowing the user to associate a diagnostic to this image.

To create a new observation, click on the `Save selected image` button in the toolbar of the *Capture Media* or *Select Media* windows or right-click on a thumbnail image or a full-size image when the multi-image mode is selected. The *Add an Exam* dialog window appears. You can then enter a diagnostic, note and change the examination date. Optional fields include diagnostic and note.

You can also create multiple observations that share the same diagnostic and related information in a single operation. In order to do so, select multiple capture images and click on the `Save selected image` button in the toolbar of the *Capture Media* or *Select Media* windows. The *Add an Exam* dialog window appears. You can then enter a diagnostic, note and change the examination date associated to the first selected image. Then an observation with the exact same diagnostic will be created for each of the remaining selected images.

The `Take a Picture and Save` button combines the acquisition of an image and the creation of an observation in a single operation.

Deleting Captured Images

Captured images are associated with the observed animal until they are explicitly deleted.

To delete captured images, select one or more images in the captured image list viewer and click on the `Delete selected image` button in the toolbar of the *Capture Media* or *Select Media* windows. The *Delete an Image* dialog window appears so you can confirm your selection before proceeding with the removal of the images.

NOTE: To disable confirmation on delete, uncheck *Always show this message* in the *Delete an Image* dialog window.

Updating an Observation

To update an observation, click on the `Update selected exam` button in the toolbar. The *Update Exam* dialog window appears allowing you to update the examination date, diagnostic and note. Each property update requires the approval of an administrator.

An entry in the audit log is added to record the change of each property.

Deleting an Observation

To delete an observation, click on the `Delete selected exam` button in the toolbar. The *Delete an Exam* dialog window appears so you can confirm your selection before proceeding with the removal of the observation. Deleting an observation requires the approval of an administrator.

An entry in the audit log is added to record the operation.

CAUTION: When deleting an observation all the depending entities are deleted as well. This includes the diagnostics and images part of this observation.

Comparing Observations

When multiple observations have been recorded, the operator may need to review the images in detail and to compare them side-by-side.

To review and compare images, click on the `Open multi-image viewer` button in the toolbar. The *Examination Viewer* dialog window appears so the operator can review and compare images side by side.

To help the operator with the comparison, selected images can be displayed in single, dual or quad mode.

Viewing Image Details

To view an image in details, click on the `Show selected image` button in the toolbar. The *Show an Image* dialog window appears allowing you to zoom in and zoom out the image.

Exporting an Image

The images can be saved into a JPEG format file. Each exported image is watermarked with specific study related information. The image tag is formatted as follows:

```
"Study Id/Animal Id/Animal Study Id/Exam Date/Day of  
Study/Laterality"
```

To save an image, click on the `Save selected image to file` button in the toolbar.

Querying the Studies Database

Keops provides a comprehensive set of in-depth searching capabilities to retrieve observations matching multi-criteria filters.

The query builder is accessed by selecting the `Search` Tab.

Building a Query

To build a query, you need to define a set of filters that match your search criteria.

To define a filter, select a property, an operator and enter a value that matches your search criteria. When entering a value in the `Value` text field, a tool-tip displays the format of the value to enter as described in the table below. Then click on the `Apply` button to add the filter to the list.

Value Type	Format
Date	'yyyy-mm-dd' Example: '2007-01-31'
Integer	[sign]digits Example: -6
Text	'[characters][%]' Example: 'Xyz123' Example with LIKE: 'X%3'

Figure 5: Property Type Formats

To combine multiple filters with logical operators (`AND`, `OR`), click on the appropriate toggle button.

The `Edit` and `Remove` buttons allow you to respectively modify or delete the selected filter from the filter list.

To execute the query, click on the `Execute query` button in the toolbar. The result can be viewed in a table format in the `Result` Tab. To view the image associated to an observation, right-click on a row in the result table and select the `View Image` menu-item. The `Messages` Tab provides detailed information about the query execution.

Exporting a Query Result

Query results can be exported into CSV format files.

To export the query result of the last query execution, click on the `Save result` button in the toolbar.

Managing the Audit Records

Keops manages modifications in compliance with FDA regulations. All changes related to study data are recorded and can be viewed in an audit log.

Recorded Operations

Most of the operations in the system are recorded for audit purpose.

Entity	Operation	Recorded
Study	Create	No
	Update	Yes
	Delete	Yes
Animal	Create	No
	Update	Yes
	Delete	Yes
Exam	Create	No
	Update	Yes
	Delete	Yes
Photo Shot	Create	No
	Delete	No
User	Create	Yes
	Update	Yes
	Delete	Yes
Lexicon Selections	Add	No
	Remove	No
Instrument	Create	Yes
	Update	Yes
	Delete	Yes
Instrument Settings	Create	No
	Update	No
	Delete	No

Figure 6: Audited Operations

Viewing the Audit Log

The audit records are listed in the `Audit` Tab. Select the `Audit` Tab to view the audit logs. The audit records can be filtered by operation type as well as by entity type.

To filter audit records, click on the `Operation` button or `Entity Type` button in the toolbar and select respectively the operations and entity type to be filtered.

Exporting the Audit Log

The audit records can be exported into a CSV format file.

To export the audit records, click on the `Save audit records` button in the toolbar.

Session Event Log

Keops records each operation executed by a Keops user during a session in an event log file. Instrument connection related events are also recorded into the log file.

The application session events are listed in the `Log` sub-tab under the `Administration Tab`.

Viewing Session Events

To view the session event records, log in as an administrator and select the `Log` sub-tab under the `Administration Tab`.

Accessing Session Events History

The `[INSTALLDIR]\logs` directory contains a history of the last 5 sessions of each connected user.

Managing the Database Server

This section gives a brief overview of the operations the administrator needs to know in order to get the Keops database server up and running.

To administer the database server, please read the *Matisse Server Administration Guide*. You can find this and other Matisse documentation on the Matisse Software web site

<http://www.matisse.com/developers/documentation/>

Creating the Keops Database

Launch Matisse Enterprise Manager and create a new database named 'keops'.

Configuring the Database

To perform these tasks with the Enterprise Manager, right click on the node that represents your host machine, and select `New Database`. Once you have entered the database name (`keops`), complete the following instructions:

1. Select the `Data Files` tab to increase the initial datafiles size to 200 Mbytes,
2. Select the `Advanced` tab to enable the database access control,
3. In the `Advanced` tab increase the server cache size to 196 Mbytes.

Starting the Database

After your database is configured, you can start it with a right-click on the database node then select `Start`. The first start will take more time as it initializes the datafiles.

Creating the Database User

You need to create the database user that Keops application is using to access the database.

Before adding a database user, you must verify that the database is online. To add a database user from the Enterprise Manager, click on the database node (`keops`) and then by clicking on the `Users` node under `Security`. Right click in the `Users` table and select `Add User`. Create a user with login name set to `keopsappli`, privileges set to `Data Read Write` and connection mode set to `Database Authentication` and enter the password of your choice.

Loading the Database Schema

The last operation consists in loading the database schema from the `KeopsApplicationSchema.odl` script file. A copy of `KeopsApplicationSchema.odl` is located in `[INSTALLDIR]\schema`

Before loading the database schema, you must verify that the database is online. To loading a database schema from the Enterprise Manager, right click on the database node (`keops`) and then select `Schema..Import ODL Schema`.

Stopping the Database Server

Before stopping a database, you must verify that there are no users connected on it. After the database has been stopped, it becomes off-line. Users are no longer able to access the data. When the database is off-line, the root or owner account can modify the database configuration file.

To stop a database from the Enterprise Manager, right click on the database node (`keops`) and then select `Stop`.

Restarting the Database Server

To restart a database from the Enterprise Manager, right click on the database node (`keops`) and then select `Start`.

NOTE: You can set the `AUTORESTART` database configuration parameter to `1` to get the database server automatically restarted after a reboot of the machine.

Backing up the Database

Matisse Database Backup tool allows users to perform full and incremental parallel backups of databases while the system in online.

There is no need to block updates during a backup, as the Matisse server keeps a snapshot of the database at the time of the beginning of the backup operation.

To backup a database from the Enterprise Manager, right click on the database node (`keops`) and then select `Tasks..Backup Database`.

Exporting and Importing Data

Keops provides scripts to export and import data in CSV files.

Exporting Data

To export data, run the export script located in `[INSTALLDIR]\data`.

Prior to running the script, you need to update the destination directory for the export of the photo shots. Edit the `PhotoShotExport.opt` file located in `[INSTALLDIR]\data` and update the destination directory of the `columnToFile` property to match your installation. The default location is `C:\Program Files\MatisseLS\Keops\data\photos`.

```
<columnToFile Directory="C:\Program Files\MatisseLS\Keops\data\photos"
      FilenameFormat="{MediaName}">Photo</columnToFile>
```

Before exporting data, you must verify that the database is online. Then open a Windows command and run `MlsKeopsExport` script with the database name as a parameter.

```
C:\> cd [INSTALLDIR]\data
C:\> MlsKeopsExport keops
Exporting Study data into Study.csv
Exporting Animal data into Animal.csv
Exporting Exam data into Exam.csv
Exporting Diagnostic data into Diagnostic.csv
Exporting PhotoShot data into PhotoShot.csv
Exporting UserAccount data into UserAccount.csv
Exporting AuditRecord data into AuditRecord.csv
Exporting Instrument data into Instrument.csv
Exporting MediaDevice data into MediaDevice.csv
Exporting InstrumentSettings data into InstrumentSettings.csv
Exporting VideoFeatureInfo data into VideoFeatureInfo.csv
Exporting UserDefinedMenuChoices data into UserDefinedMenuChoices.csv
Exporting data completed
See ExportData.log for details
C:\>
```

The CSV files are created into the `csv` directory and the images files into the `photos` directory. A log file named `ExportData.log` is produced. This file provides details about the exported data.

Importing Data

To import data, run the import script located in `[INSTALLDIR]\data`.

Before importing data, you must re-initialize the database since the database must be empty. The import script automatically loads the database schema as well as the import schema extensions. The script assumes the CSV files exist

and are located into the `csv` directory and the images files are into the `photos` directory.

To import data, open a Windows command and run `MlsKeopsImport` script with the database name as a parameter.

```
C:\> cd [INSTALLDIR]\data
C:\> MlsKeopsImport keeps
Loading Keeps Schema
Loading Animal
Loading AuditRecord
Loading Diagnostic
Loading Exam
Loading Instrument
Loading InstrumentSettings
Loading MediaDevice
Loading PhotoShot
Loading Study
Loading UserAccount
Loading UserDefinedMenuChoices
Loading VideoFeatureInfo
Linking Animal to Study
Linking Diagnostic to Exam
Linking PhotoShot to Exam
Linking Exam to Animal
Linking MediaDevice to Instrument
Linking VideoFeatureInfo to InstrumentSettings
Linking PhotoShot to Animal
Loading of data completed
See ImportData.log for details
```

A log file named `ImportData.log` is produced. This file provides details about the loaded data.

Migrating Studies

This section describes the process to migrate studies from a legacy application into Keops.

Only studies and related data are covered in this section. This includes the `Study`, `Animal`, `Exam`, `Diagnostic` and `PhotoShot` Tables.

Data Model

The following diagram describes the table format of studies to be imported into Keops.




Figure 7: Study Data Model

File Format

The CSV files must respect the following rules:

- The first row contains the column names
- Field separator is a comma (,)
- Text field delimiter is double-quotes (")
- Date Format is YYYY-MM-DD (4 digits for year, 2 digits for month, and 2 digits for day). For example January 1st, 2007 is represented as 2007-01-01
- Timestamp format is YYYY-MM-DD HH:MM:SS (4 digits for year, 2 digits for month, 2 digits for day, 2 digits for hours, 2 digits for minutes, 2 digits for seconds)

The file format for the `Study` table export reads as follows:

```
"StudyNumber", "Compound", "Species", "StartDate", "Notes"
"04200D", "PD-325000", "Rat", 2007-01-06, "this is a comment"
```

The file format for the `Animal` table export reads as follows:

```
"AnimalInitialID", "AnimalStudyID", "Group", "Sex", "Notes", "StudyN
umber_FK"
"508", "508-A", "5", "Male", "this is a comment", "04200D"
```

The file format for the `Exam` table export reads as follows:

```
"ExamID", "ExamDate", "DayOfStudy", "Notes", "StudyNumber_FK", "Anim
alInitialID_FK"
"1.JPG", 2007-01-06, 1, "this is a comment", "04200D", "508"
```

The file format for the `Diagnostic` table export reads as follows:

```
"ExamID_FK", "Laterality", "Category", "Observation"
"1.JPG", "Left Eye", "Retina", "Abnormal Vessels Nos"
```

The file format for the `PhotoShot` table export reads as follows:

```
"MediaName", "CaptureDate", "Photo", "ExamID_FK"
"1.JPG", 2002-09-16 14:21:40, 1.JPG, "1.JPG"
```

The `CaptureDate` values are timestamps. The `Photo` values represent the file name associated with the image.

Exporting Data

We assume that you are able to export from your legacy database the `study` data mapping the `Study` data model and format described above.

The export needs to produce the following CSV files:

1. `Study.csv`
2. `Animal.csv`

3. Exam.csv
4. Diagnostic.csv
5. PhotoShot.csv

In addition to the CSV files, the image files associated to the photo shots need to be produced.

Loading Data

Once you have completed the data export, copy the CSV files into Keops migration directory and copy the images files into Keops migration\photos.

To load the study data, complete the followings steps:

Re-initialize the Database

To initialize a database from the Enterprise Manager, right click on the database node (keops) and then select Re-initialize.

Load Keops Database Schema

You need to load the database schema (KeopsApplicationSchema.odl) from the Enterprise Manager. A copy of KeopsApplicationSchema.odl is located in [INSTALLDIR]\schema

Load Keops Migration Schema Extensions

You need to load the migration database schema updates (KeopsMigrationUpdates.sql) from the Enterprise Manager. A copy of KeopsMigrationUpdates.sql is located in [INSTALLDIR]\schema. To load a database schema from the Enterprise Manager, right click on the database node (keops) and then select Schema..Import DDL Schema.

Load CSV Files

The last step consists in running the migration script located in [INSTALLDIR]\migration. Open a Windows command and run MlsKeopsMigrate script with the database name as a parameter.

```
C:\> cd [INSTALLDIR]\migration
C:\> MlsKeopsMigrate keops
Loading Study
Loading Animal
Loading Exam
Loading Diagnostic
```

```
Loading PhotoShot  
Linking Diagnostic to Exam  
Linking PhotoShot to Exam  
Linking Exam to Animal  
Linking Animal to Study  
Migration completed  
See MigrateData.log for details  
C:\>
```

A log file named `MigrateData.log` is produced. This file provides details about the loaded data.

Frequently Asked Questions

This section provides quick answers to users “How do I...” questions. The FAQs are organized by category.

Configuration

How do I configure a workstation to access multiple database servers?

This is not an option of the current release. However `keopsApp.cfg` file located in Keops `config` directory contains the database server information. You can edit this file and update the database information.

Users

How do I know my user account role when logged into Keops?

The user’s role is displayed in the 4th field of the status bar at the bottom of the main window.

How do I create a Keops login without password?

All Keops users must have password.

How do I change my Keops login password?

For System users for which the username matches the Windows login name, there is no password is registered into Keops. For Application users, if you are an administrator you can change anyone password. If you are an operator, you need an administrator to log in and edit your username. You will then be able to change your password.

Acquisition Instruments

How do I know that I can capture an image using a foot pedal?

When an instrument is connected, its status indicates `External switch active` if an external switch has been configured.

User-Defined Lexicons

How do I remove the last selection in the list?

Each lexicon of closed selections must have at least one selection.

How do I enter a value that is not in the pre-defined selection list?

There is no other way around adding the value to the lexicon first.

Studies

How do I compare side-by-side observations from 2 different animals?

This is not an option of the current release.

How do I print an observation including image and diagnostic?

There is no print function the current release. However you can export an observation image which is tagged with a watermark detailing the observation origin. You can also export the result of a search into an Excel type file.

Audit Records

How do I remove a record from the audit log?

You cannot remove any audit record for the audit log.

How do I reset the audit log?

The audit log cannot be reset or cleared of any record.